



STAND-UP BAGS

7 SIZES AVAILABLE



OPENING PULL TABS

- for easy opening
- non-slip textured Tyvek®

CONVENIENT WRITE-ON STRIP

- makes labeling quick and easy
- clear sample identification
- for organization & traceability

GUSSETED BAG BOTTOMS FOR SELF-STANDING

- no need for expensive and bulky bottles
- allows user to perform functions with both hands

STAND-UP BAG WHIRL-PAK® STAND-UP



TRUSTED, PUNCTURE-PROOF CLOSURE TABS

- tape extends beyond wire to eliminate sharp ends and provide a sturdy grip
- wire tabs hold the top open for easy access to add to or remove sample
- worry-free, no-leak closure for sample integrity

+1

POST-MANUFACTURING STERILIZATION

- increased sample integrity
- perforated closure to maintain sterility



BLENDED, VIRGIN POLYETHYLENE FILM

- exceptional clarity
- superior strength
- sterilized after manufacturing
- FDA approved for contact

more efficient than bottles



smarter sustainability

LESS STORAGE SPACE • EASY TO SHIP • USER FRIENDLY • COST EFFECTIVE



+ ideal for use in water, environmental, protein and more.

Whirl-Pak® sterilized sampling bags are R-Nase, D-Nase, Pyrogen and BPH free and comply with regulatory and industrial requirements:





500 BAGS VS. 500 BOTTLES


GUARANTEED STERILE
Manufactured according to
quality management systems
ISO9001:2015

- o Uses less shelf space
- o Easier disposal
- o Durable
- o Lower shipping costs
- o Built-in closure
- o Lower total cost



SOLIDS OR LIQUIDS - LEAK-PROOF CLOSURE

STEP 1




Label your sample if using a write-on bag.

STEP 2



To open, tear off the top along the perforation.

STEP 3



Use the Tyvek® pull tabs to open without contamination.

STEP 4




Place sample in bag. Perfect for liquids or solids.

STEP 5A




Hold bag by the tabs and whirl 3 to 4 times to create a seal.

STEP 5B



Or, fold the tab over tightly 3 to 4 times to close.

STEP 6A



Secure closure, bend wire ends for solids.

STEP 6B



Or, bring the two wire ends together to twist, tie and tuck!

TO LEARN MORE **CONTACT**



WHIRL-PAK[®]
Filtration Group



whirl-pak.com



+1.920.568.5616



contact@whirl-pak.com